



Quality First
Develop Harmoniously

EN ISO 13485:2016

CE

FDA



CITEC™ CLEANING BRUSH FOR ENDOSCOPE


Date of issue: 2024-08-12

Cleaning Brush for Endoscope

Features

- Complete product range for endoscope channel cleaning and channel opening cleaning.
- Highly resistant nylon made bristles, ensure maximum clean effect.
- Flexible and kink resistant sheath allows smooth insertion throughout the endoscope channel.
- Round tip brush end prevents from accidental damage or scratches on endoscope.
- Single-use device reduces the risk of cross-contamination.
- Ergonomic handle design for easy manipulation.



Cleaning Brush for Endoscope Channel



Ref. Code	L(mm)	D(mm)	L1(mm)	D1(mm)	L2(mm)	D2(mm)	Channel Size(mm)	Units/ Box
CIT-CB-0303-2300-C	2300	1.7	20	3	20	3	2.0 ~ 2.8	100
CIT-CB-0505-2300-C	2300	1.7	20	5	20	5	≥2.0	100
CIT-CB-0606-2300-C	2300	1.7	20	6	20	6	≥2.0	100



Remark: Double bristle ends.






Ref. Code	L(mm)	D(mm)	L1(mm)	D1(mm)	D2(mm)	L2(mm)	D3(mm)	D4(mm)	Channel Size(mm)	Units/ Box
CIT-CB-0507-2300-C	2300	1.7	20	5	7	20	7	5	≥2.0	100

Remark: Double tapered bristle ends.



Ref. Code	L(mm)	D(mm)	L1(mm)	D1(mm)	L2(mm)	D2(mm)	Channel Size(mm)	Units/ Box
CIT-CB-0505-2400-C	2400	1.7	20	5	125	5	≥2.0	100


Remark : Double ends,one bristle end and the other pipe end.




Ref. Code	L(mm)	D(mm)	L1(mm)	D1(mm)	Channel Size(mm)	Units/ Box
CIT-CB-2.4-900-C	900	0.7	5	2.4	1.0-1.5	100
CIT-CB-2.4-1650-C	1650	0.7	5	2.4	1.0-1.5	100

Remark: Single bristle end, with flexible metal sheath.


Cleaning Brush for Channel Opening



Ref. Code	L(mm)	L1(mm)	D1(mm)	Units/ Box
CIT-CB-V1	55	35	10	200




Ref. Code	L(mm)	L1(mm)	D1(mm)	L2(mm)	D2(mm)	Units/ Box
CIT-CB-V2	155	20	5	30	10	200
CIT-CB-V4	155	20	5	30	11	200





Ref. Code	L(mm)	L1(mm)	D1(mm)	L2(mm)	D2(mm)	Units/ Box
CIT-CB-V3	190	40	5	40	12	200

Combination Cleaning Brush for Endoscope



Ref. Code	L(mm)	D(mm)	L1(mm)	D1(mm)	L2(mm)	D2(mm)	Channel Size(mm)	Units/ Box
CIT-CB-0510-2300-C	2300	1.7	20	5	45	10	≥2.0	100


Remark: Double bristle ends. One end for endoscope channel cleaning and the other end for channel opening cleaning.




Ref. Code	L(mm)	D(mm)	L1(mm)	D1(mm)	L2(mm)	D2(mm)	Channel Size(mm)	Units/ Box
CIT-CB-0510-2390-C	2390	1.7	20	5	35	10	≥2.0	100

Remark: Double bristle ends. One end for endoscope channel cleaning and the other end for channel opening cleaning.

Cleaning Brush Kit




CIT-CB-KIT1




CIT-CB-KIT2



CIT-CB-KIT3



CIT-CB-KIT4



CIT-CB-KIT5

Ref. Code	Composition	Kits/ Box
CIT-CB-KIT1	Cleaning Brush for Channel Opening (CIT-CB-V1)	100
	Cleaning Brush for Endoscope Channel (CIT-CB-0505-2300-C)	
CIT-CB-KIT2	Cleaning Brush for Channel Opening (CIT-CB-V2)	100
	Cleaning Brush for Endoscope Channel (CIT-CB-0505-2300-C)	
CIT-CB-KIT3	Cleaning Brush for Channel Opening (CIT-CB-V3)	100
	Cleaning Brush for Endoscope Channel (CIT-CB-0507-2300-C)	
CIT-CB-KIT4	Cleaning Brush for Channel Opening (CIT-CB-V4)	100
	Cleaning Brush for Endoscope Channel (CIT-CB-0606-2300-C)	
CIT-CB-KIT5	Cleaning Brush for Channel Opening (CIT-CB-V4)	100
	Cleaning Brush for Endoscope Channel (CIT-CB-0505-2400-C)	

Remark: We can provide the customized kit with different composition.

