


Disposable Esophageal Bougie



Intended use

- Used for dilatation of esophageal stenosis.

Features

- Made of medical macromolecule, with proper rigidity.
- Equipped with 2 radiopaque marks on each bougie, easy to observe under X-ray.
- Various bougie sizes, meet different requirements.

Product code	Bougie diameter (mm)	Units/Box
CIT-EB-05	5	20
CIT-EB-06	6	20
CIT-EB-07	7	20
CIT-EB-08	8	20
CIT-EB-09	9	20
CIT-EB-10	10	15
CIT-EB-11	11	15
CIT-EB-12	12	15

Product code	Bougie diameter (mm)	Units/Box
CIT-EB-13	13	15
CIT-EB-14	14	15
CIT-EB-15	15	10
CIT-EB-16	16	10
CIT-EB-17	17	10
CIT-EB-18	18	10
CIT-EB-19	19	10
CIT-EB-20	20	10