


Disposable Veress Needle



Features:

- + Retractable style safety mechanism, spring loaded tip protector.
- + Visible needle position with red indicator.
- + Sharp tip for easy abdominal puncture.
- + Ergonomic wing or round handle provides comfortable and convenient manipulation.
- + 120/150 mm length available for different procedure requirements.
- + High quality and cost effective.

| Code No. | Specifications | Units/ Box |
|---------------|---------------------|------------|
| CIT-QFZ-120 | Wing handle, 120mm | 50 |
| CIT-QFZ-150 | Wing handle, 150mm | 50 |
| CIT-QFZ-120-1 | Round handle, 120mm | 50 |
| CIT-QFZ-150-1 | Round handle, 150mm | 50 |