


# Endoscope Channels Flushing Tubing Set



## Intended Use

The tubing set is used for channel flushing of flexible immersible endoscopes during manual cleaning process.

## Features

- Single use tubing set, reduce the working burden of cleaning.
- With three individual clamps, realize easy control the channel cleaning.
- If necessary, pump tube can be cut into two parts and reconnect with the connector.

Ref. Code	Description	Units/Box
CIT-FTS-1	Connector *1, Block Adapter*1, Water Jet Connector*1, Length of pump tube:1.2m, for Olympus Endoscopes	15
CIT-FTS-2	Connector *1, Block Adapter*1, Water Jet Connector*1, Length of pump tube:2.0m, for Olympus Endoscopes	15

