

Disposable Bipolar Instruments



Features:

- + Two kinds of Jaws configuration available, dual action or single action.
- + Jaws rotate 360° freely via the rotation knob.
- + Ergonomic handle design to provide comfort gripping and manipulation.
- + Bigger jaws opening for both dissector and grasper.

Code No.	Product	Spec.(mm)	Units/ Box
07F-004S	Maryland, single Action	5 x 330	5
07F-004D	Maryland, double Action	5 x 330	5
07F-005S	Endo grasper, single Action	5 x 330	5
07F-005D	Endo grasper, double Action	5 x 330	5

Different Jaws:



Maryland, single action



Maryland, double action



Endo grasper, single action



Endo grasper, double action